

Instructions for Whitening with Style White



1
Brush teeth prior to wearing the tray. Make sure the trays are clean and dry before applying the gel. Place a small dab of tooth whitening gel half way up the front of each tooth in the tray. We recommend that you only whiten the teeth which are visible when you smile. Do not spread the gel around. Cap and retain syringe for next treatment until empty.



2
Insert whitening tray in the mouth over the teeth. Seat the tray firmly against the teeth.



3
Gently remove excess gel that may flow beyond the tray's edge on to the gums. To do this, hold the tray up with one hand and wipe the excess gel off the gums with your finger or a tissue. Wear the tray with the gel for 30 minutes once a day.

Do Not wear overnight.



4
Remove the tray and rinse off your teeth. Brush away any remaining gel. Clean tray with toothbrush and cold water. Store tray and gel away from heat and direct sunlight.

**Have more questions?
View Frequently Asked Questions at
www.stylewhitesmiles.com**

**Style
White**
INSTRUCTIONS

“Keeping the smile of a lifetime.”

One of the biggest issues with teeth whitening is the tendency for teeth to turn dark again over time. This is true for all whitening systems. The indulgences of coffee, tea, red wine, tobacco, berries and many more products are the very things that will darken your youthful smile over time. In order to maintain your beautiful white smile you will have to whiten your teeth every 6-24 months depending on your habits. This is accomplished in exactly the same way you initially whitened them. You may purchase additional whitening gel from Style White. Visit us at www.stylewhitesmiles.com

IF YOUR WHITENING GEL IS OVER 2 YEARS OLD IT HAS EXPIRED AND IS NO LONGER EFFECTIVE.